




Stop the Insanity

CS-023
Designed by
Steven Bressan



Drawing is Not to Scale!

Start Position:

String 1: Standing in Box A, facing downrange, hands relaxed at sides.

String 2: Standing in Box A, facing uprange, hands above respective shoulders.

Procedure:

String 1: At start signal engage T1 thru T3 with only two (2) rounds each in any order, perform a mandatory reload and re-engage T1 thru T3 with only two (2) rounds each in any order.

String 2: At start signal, assume mandatory prone position and engage T1 thru T3 with only two (2) rounds each in any order.

Scoring: Shots Limited

Targets: 3 NRA D-1 paper

Rounds: 18 maximum

Start: Audible

Stop: Last shot

Penalties: Per the ICORE rulebook. Foot faults, extra shots, extra hits, procedurals and misses +5 seconds per occurrence. Failure to go prone incurs a single procedural.

Stage Setup: Place targets per drawing dimensions. Set targets T1 thru T3 at standard height with top of target 1.65 m (5 feet 6 inch) +/- 5cm (2 inches) from ground. Alignment is shown to centerlines of boxes and targets, unless otherwise noted.

Shooting boxes are approximately 1m x 1m (3 feet x 3 feet).

CS-023

Shots Limited

Stop The Insanity



R.O. _____

Target	A	B	C	M	Hits
T1					6
T2					6
T3					6

Shooter _____

String 1 Time

String 2 Time

Total Hits: + + + = 18

Time Factor: x0 x1 x2 x5

+ + =

Extra Shots X5=

Extra Hits X5=

Procedurals X5=

Notes

Subtotal Time

Hit Time Added

+ Extra Shots

+ Extra Hits

+ Procedurals

Total

Name _____

ICORE Number _____

Date _____ Club _____

Limited / Open / Retro (Circle one)